

Deluxe Anesthesia Cart #6456



Specifications

Construction:

- Steel

Features include:

- Raised back rail system
- 6-Compartment tilt bin organizer
- Tape and label dispenser
- Large utility basket
- Catheter box
- Large utility basket
- 6-Socket hospital grade outlet with 15-foot cord and cord wrap
- Plastic waste container
- Large sharps container and bracket
- 3-inch Exchange tray (2)
- 6-inch Exchange tray (2)
- 12-inch Adjustable drawer divider set
- Tall cabinet with key locking system
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Plastic removable top
- Aluminum push handles (4)
- Ball-bearing drawer slides
- Double-walled drawer faces with no rough edges
- Cart is delivered wheels assembled
- Certified "Latex Free"
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

Drawer configuration:

- 4 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)
- 1x12-inch (30.5 cm)

Drawer dimensions:

- 23" W x 17" D (58.4 cm x 43.2 cm)

Twelve available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Assembled cart dimensions:

- 41.5" H x 32" W x 22" D (105 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 171 pounds (77.8 kg)
- Approx. shipping weight: 206 pounds (93.7 kg)