

Three Drawer Infection Control Cart with Breakaway Lock Specialty Package #6501

Specifications



Construction:

- All steel

Features include:

- Coat and hat rack
- Utility hooks, set-4
- Stainless steel top rail/push handle
- Tall cabinet with stainless steel breakaway locking system
- Breakaway lock seals, pkg-50
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Removable plastic top
- Ball-bearing drawer slides
- Double-walled drawer faces with no rough edges
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

Twelve available color choices:

- Yellow (-YL), Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Drawer configuration:

- 1 x 6-inch (15.2 cm)
- 2 x 12-inch (30.5 cm)

Assembled cart dimensions:

- 44.5" H x 32" W x 22" D (113 cm H x 81.3 cm W x 55.9 cm D)
- With coat and hat rack 62.5" (158.8 cm)
- Cart weight: 145 pounds (66 kg)
- Shipping weight: 180 pounds (81.9 kg)