



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

OptimAL Line, Aluminum Six Drawer MRI-Compatible Emergency Cart, Breakaway Lock, Standard Package #MR6B



Specifications

Construction:

- Non-ferrous painted aluminum cabinet and drawers

Features:

- Designed for compatibility in MRI suites
- Tested for safe use in MRI suites by independent testing authority (copies of report available)
- Cabinet painted gloss white with MRI identification sticker applied
- 13 drawer colors available (see below)
- MRI-compatible accessories available
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel – one with brake, one directional
- Stainless steel breakaway locking system
- Breakaway lock seals, numbered, yellow, pkg-50
- Full wrap-around vinyl bumper
- Stainless steel pull-out shelf
- Replaceable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- Certified "Latex Free"
- Five year warranty
- Hard-baked powder coat finish

Thirteen Available Color Choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Grey (-TXGY), Hammertone Green (-TXGN) and Sand (-SD).

Assembled Cart Dimensions

- 45.75" H x 30.5" W x 23.75" D
(116.2 cm H x 77.5 cm W x 60.3 cm D)
- Cart net weight: 126 pounds (57.3 kg)
- Approx. shipping weight: 146 pounds (66.4 kg)

Drawer Configuration:

- 1 x 4" (10.2 cm)
- 2 x 3.25" (8.3 cm)
- 2 x 6.5" (16.5 cm)
- 1 x 9.75" (24.8 cm)

Internal Drawer Dimensions:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)