

Multi-Treatment Cart Standard Package #6200

Specifications



Construction:

- All steel

Features include:

- Tall Cabinet (with key locking system for drawers and doors)
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) premium single-wheel caster set, full swivel – one with brake, one directional
- Ball-Bearing Drawer Slides
- Double-walled drawer faces with no rough edges
- Removable Plastic Top
- Aluminum Push Handle
- Storage Compartment: 21-inch (53.3 cm) high opening, with adjustable internal shelf
- Cart is delivered fully assembled
- Certified "Latex Free"
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Drawer configuration:

- 1 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)

Assembled Cart Dimensions:

- 41.5" H x 32" W x 22" D (105.4 cm H x 81.3 cm W x 55.9 cm D)
- Cart Net Weight: 115 pounds (52.2 kg)
- Approx. Shipping Weight: 150 pounds (68.0 kg)