

MED-BIN CARTS

SL36BIN5 (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 36 - 5" bins, Standard Package

SL28BIN5 (not shown)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 28 - 5" bins, Standard Package

SL54BIN3 (pictured)
Standard Line Med-Bin Cart with a BEST® lock on cabinet, pull-out shelf. Capacity of 54 - 3.5" bins, Standard Package

SL42BIN3 (not shown)
Standard Line Med-Bin Cart with a BEST® lock on cabinet, pull-out shelf. Capacity of 42 - 3.5" bins, Standard Package



Bin dividers & labels included with each cart

SL36BIN5XL (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 36 extra long 5.6" bins, Standard Package

SL24BIN5XL (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 24 extra long 5.6" bins, Standard Package

SL18BIN5XL (pictured)
Standard Line Med-Bin Cart with BEST® lock on cabinet, pull-out shelf. Includes 18 extra long 5.6" bins, Standard Package



Bin transfer cart options also available



Bin drawer with XL bins

Bin dividers & labels included with each cart



40442
Replacement
3.5" Bins
2.75"H x 3.5"W x 9.5"D

680533
Bin Label for
3.5" Bins
(150 Labels)

680513
10 Bin Dividers
for 3.5" Bins



40446
Replacement
5" Bins
2.75"H x 5"W x 9.5"D

680535
Bin Labels for
5" Bins
(100 Labels)

680515
10 Bin Dividers
for 5" Bins



40433
Replacement
5.6" Extra Long Bins
2.75"H x 5.6"W x 18"D

680537
Bin Labels for
5.6" Bins
(210 Labels)

40434
20 Bin Dividers
for 5.6" Bins

VL36BIN3 (pictured)
Value Line Med-Bin Cart with
key lock on cabinet, internal
narcotics box, pull-out shelf.
Includes 36 - 3.5" bins,
Standard Package




**Mechanical keyless entry
available on
Value Line Carts**



VL24BIN5 (pictured)
Value Line Med-Bin Cart with key lock on
cabinet, internal narcotics box, pull-out shelf.
Includes 24 - 5" bins, Standard Package

VL24MBIN5 (not shown)
Value Line Med-Bin Cart with keyless lock on
cabinet, internal narcotics box, pull-out shelf.
Includes 24 - 5" bins, Standard Package

VL36MBIN3 (not shown)
Value Line Med-Bin Cart with keyless lock on
cabinet, internal narcotics box, pull-out shelf.
Includes 36 - 3.5" bins, Standard Package



Bin Drawer shown with 3.5" bins

Harloff standard bin carts are available
with 5" and 3.5" bins. They extend
8.5" into the drawer and have a
small tab that fits under a
metal edge. There is 7" of
open space behind



VL12BIN5 (pictured)
Value Line Med-Bin Cart with
key lock on cabinet, internal
narcotics box, pull-out shelf. Includes
12 - 5" bins, Standard Package
(shown with a Specialty Package)

VL18BIN3 (pictured)
Value Line Med-Bin Cart with
key lock on cabinet, internal
narcotics box, pull-out shelf.
Includes 18 - 3.5" bins,
Standard Package

**Accessory packages available
for all cart styles**