



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Medical Supply Visual Bin Storage Cart, with Louvered Panel for Akro-Mils Storage Bins; Hinged Plexiglas Door, Key Lock MS4372LVR



*Bins not included.

Thirteen available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY) Hammertone Green (-HTGN) and Sand (-SD).

Assembled Cart Dimensions

- 71.5" H x 43" W x 26.5" D (181.6 cm H x 58.4 cm W x 67.3 cm D)
- Cart weight: 258 pounds (117.3 kg)
- Shipping weight: 370 pounds (167.8 kg)

Advantages and Benefits

Harloff's MS4372LVR Visual Bin Storage Cart is designed for use where space for storage in the OR suite and SPD sterile storage area is at a premium and nurses need to locate supplies in a hurry.

The cart provides flexible storage and transport of sterile medical supplies while helping facilities secure and manage inventory to cut costs, decrease waste, improve productivity and enhance patient care.

It's unique feature is a louvered panel inside the locking cabinet which holds up to 42 Akro-Mills AkroBins® (bins not included). Use the multiple sizes and colors of these bins to organize supplies in standardized locations that can be delivered to wherever they are needed, increasing workspace efficiency and improving workflow.

Quick View Plexiglas panels allow staff to quickly check stock.

Cart is delivered fully assembled and ready to use.

Make your OR/SPD supply chain run smoothly with Harloff's Medical Supply Visual Storage Cart.

Specifications

Construction:

- Painted 20-gauge powder-coated steel cabinet, Plexiglas, steel louvered panel

Superior Features:

- Enclosed and locking cabinet keeps items secure.
- Standard clear front doors allow visual inspection of supplies.
- Built using powder-coated steel for durability and easy to clean surfaces

Features:

- Louvered Back Panel holds multiple sizes of AKROBINS® (bins not included)
- Clear front key locking doors with Plexiglas inserts
- Key locking doors open flush to sides of cabinet for unobstructed access
- Integrated side push handles
- Full wrap-around vinyl bumper
- 5" (12.7 cm) premium rubber full-swivel ball-bearing casters (1 with brake, 1 tracking)
- Cart is delivered with wheels installed
- Latex free
- Twelve Year Warranty
- Hard-baked powder coat finish