



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Six Drawer Mini Line Emergency Cart with Wrap-around Bumper and 5” Wheels, with Break-away Lock and Specialty Package #3256B-EMG



Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Includes MINI-BPRWHL- a wrap-around stabilizing bumper and 5-inch (12.7 cm) premium caster set, full swivel – two with brakes
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Break-away locking system
- 3” (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-Prong IV pole (extends to 68.75” H) (680401)
- Oxygen tank holder (680403)
- Cardiac board and brackets (680406)
- Aluminum push handle / side rail (68161)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Sixteen available color choices:

***Antimicrobial Finish Available**

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 15” W x 16.75” D (38.1cm x 42.5cm)

Assembled Cart Dimensions

- 40” H x 25” W x 22”D (101.6 cm H x 63.5 cm W x 55.9 cm D)
- With accessories: 68.75” H x 27” W x 22” D (174.6 cm H x 68.5 cm W x 55.9 cm D)
- Cart net weight: 118 pounds (53.5 kg)