



HARLOFF
650 Ford Street
Colorado Springs, CO 80915

How to specify and order UNIVERSAL LINE CARTS

Universal Line carts are identified by a standardized sequence based on Body Size, Locking Type, Number of Drawers, additional package, and any other special designations.

Body Size/Designation (Working height)	Locking Type	Drawers	Package
810 – Super (45")	B=Breakaway	7 = Standard Seven for Body Size	EMG Emergency
809 – Tall (41 ½")	K=Keylock	6= Standard Six for Body Size	ANS Anesthesia
808 – Medium(38")	M=Mechanical Combo	5= Standard Five for Body Size	ANS2 Deluxe Anesthesia
		4= Standard Four for Body Size	INF Infection Control
		3= Standard Three for Body Size (Tall and Medium bodies only)	

For example, a 809B6-EMG is a tall cart with a breakaway lock, the standard six drawer configuration for the tall cart (4-3.25", 1-6.5", and 1-9.75"), and the EMG emergency package included.

Standard Drawer Specifications – Selected Models

Base Product number	Standard Drawer Configuration
810B7/810K7/810M7	5-3.25" (8.3 cm), 1-6.5" (16.5 cm), 1-9.75" (24.8 cm)
809B6/809K6/809M6	4-3.25" (8.3 cm), 1-6.5" (16.5 cm), 1-9.75" (24.8 cm)
809K4	1-3.25" (8.3 cm), 1-6.5" (16.5 cm), 2-9.75" (24.8 cm)
809K3	3-9.75" (24.8 cm)
808B5/808K5/808M5	3-3.25" (8.3 cm), 1-6.5" (16.5 cm), 1-9.75" (24.8 cm)
808K4	1-3.25" (8.3 cm), 2-6.5" (16.5 cm), 1-9.75" (24.8 cm)
808K3	1-6.5" (16.5 cm), 2-9.75" (24.8 cm)

To specify a custom drawer configuration, Harloff uses the following convention:

Place the letter "C" after the body, lock and drawer code. For example, 808K6C identifies a Medium Keylocking cart that requires six custom configured drawers within the 8 available 3.25" "slots". Drawers must be specified by customer in sequence top to bottom. For example, this customer may specify 1-6.5" top drawer, 4-3.25" middle drawers, and 1-6.5" bottom drawer to fill the available drawer space.

Standard Overall Dimensions Without Accessories

Body Style	Dimensions (Height, Width and depth)	Estimated Net Weight
810 Carts/7 drawer	48" H x 31" W x 24.75" D (122 cm H x 78.7 W cm x 63 cm D)	162 lbs/ 74 kg
809 carts/6 Drawer	44.5" H x 31" W x 24.75" D (113 cm H x 78.7 W cm x 63 cm D)	148 lbs/ 67.3 kg
808 carts/5 Drawer	41" H x 31" W x 24.75" D (104.1 cm H x 78.7 W cm x 63 cm D)	135 lbs/61.4 kg



719.637.0300

medicalsales@harloff.com

www.harloff.com



**UNIVERSAL LINE
MEDICAL CARTS FOR
CLINICAL ENVIRONMENTS**



THE MEDICAL CART REDEFINED

- Modern Design at an Affordable Price
- Built for the Full Range of Clinical Department Needs
- Accommodates Our Full Range of Accessories and Drawer Dividers
- Standard and Custom Drawer Configurations Available
- Multiple Security Systems for All Requirements
- Proudly Made in the USA



Storage, Delivery and Security Solutions since 1951

www.harloff.com

719-637-0300

medicalsales@harloff.com

EMERGENCY CARTS

Six standard models in three cabinet heights with a wide variety of drawer configurations. Base carts include the following features:

- Breakaway lock with 100 plastic numbered seals • Internal pull-out shelf • Top rail and push handle • Plastic top
- Five inch full swivel casters, one with brake, one directional • Protective extended pontoon bumpers for stability
- Ball bearing drawer slides • Pre-specified drawer configuration

EMG emergency package upgrade includes: D/E oxygen cylinder holder

- IV pole • Rear-mounted cardiac board • Utility hooks (2)
- Three inch full drawer trays with adjustable dividers (2)



808B5
Medium Five Drawer



809B6-EMG
Tall Six Drawer with EMG Package



810B7
Super Seven Drawer

Treatment and Infection Control Carts

The Universal Line offers three cabinet heights, selection of locking systems, and a wide variety of standard drawer configurations to meet your department's storage and security requirements. A sample of these carts is shown below. Call Harloff at **800-433-4064** or e-mail to medicalsales@harloff.com for a complete line sheet.



808K3-INF
Medium Three Drawer
with INF Package



808K5
Medium Five Drawer



809M6
Tall Six Drawer



808K3
Medium Three Drawer

ANESTHESIA CARTS

Six standard models in three cabinet heights and multiple standard accessory packages. Base carts include the following features:

- Selection of key lock or mechanical combination lock • Internal pull-out shelf • Top rail and push handle • Plastic top
- Five inch full swivel casters, one with brake, one directional • Protective extended pontoon bumpers for stability
- Ball bearing drawer slides • Pre-specified drawer configuration



808M5
Medium Five Drawer



809K6
Tall Six Drawer



810M7
Super Seven Drawer



808M5-ANS
Medium Five Drawer
with ANS Package



809M6-ANS
Tall Six Drawer
with ANS Package



810K7-ANS
Super Seven Drawer
with ANS Package

ANS anesthesia package upgrade include: Three tier back rail system • Four side mount rails • Six socket hospital grade electrical outlet with 15 foot cord and cord wrap.

ANS2 deluxe anesthesia package upgrade includes all ANS package items plus six compartment tilt bin • Label holder • Catheter box • Three gallon sharps container • Three gallon waste container • Three inch full drawer tray with dividers • Six inch exchange trays (2), Nine inch drawer divider set

All carts shown with standard drawer configurations, custom configurations also available.