

V-SERIES MEDICAL CARTS
FOR CLINICAL ENVIRONMENTS



V30-2K



V24-3K



V30-6EMG



VALUE PRICED MEDICAL CARTS
DESIGNED FOR HEALTHCARE FACILITIES
UNDER BUDGET PRESSURE.



V30-6K

HARLOFF
www.harloff.com

Quality Construction at a Value Price

EMERGENCY CARTS

Perfect for Critical Supplies
used in Emergency Situations



INFECTION CONTROL/ISOLATION CARTS

Large Drawers Accommodate PPE Supplies



TREATMENT AND PROCEDURE CARTS

Ideal for use in a Variety of
Clinical Applications



V-SERIES CARTS – COMPETITIVELY PRICED WITH THE SAME HARLOFF QUALITY

- Extended wheel base with projected front casters to prevent tipping
- Seamless painted steel top work surface
- Tall or short cabinet heights available
- Carts available with break-away lock bar or key lock
- Integrated drawer pull handles
- Full-extension ball-bearing drawer slides
- Five inch full swivel casters, two with brake
- Ships with wheels installed and ready to use
- Twelve year manufacturer's warranty
- Latex free
- Durable, hard-baked powder coat finish



Call Harloff at **800-433-4064** or e-mail to medicalsales@harloff.com for a complete line sheet.

719.637.0300

medicalsales@harloff.com

www.harloff.com

V-SERIES ACCESSORIES

V-Series carts are pre-drilled to accommodate V-Series specific accessories

V-Series accessories are customer installed



VS-TOPRAIL Top Rail and Push Handle, Painted Steel



VS-DRPSHELF Side Mounted Drop Shelf without Lip



VS-CARDBRD Cardiac Board with 3 Brackets



VS-AL02 Oxygen Tank Holder, for D/E Cylinders



VS-IV2 IV Pole, 2 Prong



68530-P1 3" Drawer Divider Tray with Adjustable Dividers



681502 6" Exchange Tray with Adjustable Dividers, Fits 2 per drawer



680503 9" and 12" Drawer Divider Set with 3 Long and 2 Short Dividers



VS-3GWASTE Plastic Waste Container Mounting Bracket without Cover, 3 Gal.



VS-8QWASTE Plastic Waste Container without Cover, 2 gallon



VS-3GSHARPS 3 Gal Bemis® Sharps Container with Wire Holder and Mounting Bracket



VS-UNIVSHARP Sharps Bracket, Universal



VS-GLVBOX Glove Box Holder Plastic for 684810/40835 Sharps Holder



VS-PLTOP Top, Plastic (upgrade for factory installation)

| | MODEL | DRAWER CONFIGURATION | OVERALL DIMENSIONS | WORK SURFACE HEIGHT | WEIGHT |
|--------------|-----------------|--|--|---------------------|-------------------|
| Short Height | V24-3K | 1-6" (15.2 cm), 2-9" (22.9 cm) | 34" H x 29.5" W x 22" D (86.4 cm H x 74.9 cm W x 55.9 cm D) | 34" | 88 lbs (39.9 kg) |
| | V24-4K | 2-3" (7.6 cm), 1-6" (15.2 cm), 1-12" (30.5 cm) | 34" H x 29.5" W x 22" D (86.4 cm H x 74.9 cm W x 55.9 cm D) | 34" | 98 lbs (44.5 kg) |
| | V24-4B | 2-3" (7.6 cm), 1-6" (15.2 cm), 1-12" (30.5 cm) | 34" H x 29.5" W x 22" D (86.4 cm H x 74.9 cm W x 55.9 cm D) | 34" | 98 lbs (44.5 kg) |
| | V24-4EMG | 2-3" (7.6 cm), 1-6" (15.2 cm), 1-12" (30.5 cm) | 61.5" H x 33.25" W x 22.5" D (156.2 cm H x 96.5 cm W x 55.9 cm D) height includes IV pole - extends to 79" H (200.6 cm) | 34" | 110 lbs (49.9 kg) |
| | V24-5K | 3-3" (7.6 cm), 1-6" (15.2 cm), 1-9" (22.9 cm) | 34" H x 29.5" W x 22" D (86.4 cm H x 74.9 cm W x 55.9 cm D) | 34" | 101 lbs (45.8 kg) |
| Tall Height | V30-2K | 1-3" (7.6 cm), 1-6" (15.2 cm) Drawers and 1-21" (53.34 cm) Storage Compartment | 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D) | 40" | 90 lbs (40.8 kg) |
| | V30-3K | 1-6" (15.2 cm), 2-12" (30.5 cm) | 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D) | 40" | 102 lbs (46.3 kg) |
| | V30-4K | 1-3" (7.6 cm), 3-9" (22.9 cm) | 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D) | 40" | 105 lbs (47.6 kg) |
| | V30-5K | 2-3" (7.6 cm), 1-6" (15.2 cm), 2-9" (22.9 cm) | 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D) | 40" | 115 lbs (52.2 kg) |
| | V30-6K | 4-3" (7.6 cm), 1-6" (15.2 cm), 1-12" (30.5 cm) | 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D) | 40" | 118 lbs (53.5 kg) |
| | V30-6B | 4-3" (7.6 cm), 1-6" (15.2 cm), 1-12" (30.5 cm) | 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D) | 40" | 118 lbs (53.5 kg) |
| | V30-6EMG | 4-3" (7.6 cm), 1-6" (15.2 cm), 1-12" (30.5 cm) | 67" H x 33.25" W x 22.5" D (170 cm H x 84.5 cm W x 57.2 cm D) height includes IV Pole - extends to 85" H (216 cm) | 40" | 138 lbs (62.6 kg) |