

WOOD LOOK CARTS & CABINETS

Skirting Conceals Wheels

Locking Narcotics Storage

In Wall Medication Cabinet

HARLOFF
www.harloff.com

WOOD LAMINATE CARTS

Harloff's line of wood look carts are designed to support Culture Change initiatives. Even though these carts and cabinets look like furniture, the internal structure is engineered to meet industry standards for security and function.

Features Include:

- All aluminum internal structure
- Proven locking system
- Internal narcotics box
- Pushable
- Skirting conceals 4" ball-bearing wheels
- MDF panels covered with a special heat pressed thermal foil laminate
- Choice of laminate finishes

Designed to accommodate a variety of medication packaging.



WL294BOX-CM

WL480PC-CM



WL6K-CM



WL240PCN-SO



WL108BOX-SO

IN WALL CABINET-SINGLE UNIT



WL2780TRMPKG



CHOICE OF CABINET DOOR STYLES

IN WALL CABINET-DOUBLE UNIT



WL2780DBREC-CM

Residents can share a cabinet or have one of their own. Our double unit adds room for the storage of treatment supplies in addition to medications.

IN-ROOM RESIDENT MEDICATION CABINETS

Functional solution for resident medication storage

- Options for 1-2 residents
- Storage for medications, liquids and treatment supplies
- Painted steel cabinets
- Door is MDF panel covered with a special heat pressed thermal foil laminate
- One hand open and close
- Predrilled holes for mounting
- Mounting hardware included
- Choice of laminate finishes

WALL MOUNT CABINET



WL2790-CM

Cabinet accommodates treatment supplies, liquids and medications. Narcotics storage option is available.



CABINET DOOR STYLES

WL2790 TRIM OPTION



CABINET FRAME

CABINET WRAP

WALL MOUNT CABINET



WL2717

The WL2717 is perfect for rooms with wardrobe closets.



650 Ford Street
Colorado Springs, CO 80915



WL240PCN-S0
Punch Card Cabinet

Harloff manufactures a complete line of wood furniture for use in Long Term Care Facilities. These products support culture change initiatives. The Harloff product line features emergency, treatment, medication carts, and in-room resident medication cabinets.



WV6400-CM
Emergency Cart



WV6450-OAK
Treatment Cart



WV540PC-CM
Medication Cart



**Rear Mounted
Accessories Preserve
the Look of Furniture**

Call a member of our Long Term Care Team at **719.637.0300 ext. 108** or **129** or email us at medicalsales@harloff.com.

