



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Preconfigured Double Column Medical Storage Cart with Tambour Door MS8140-VASC

The MS8140-VASC is a preconfigured double column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8140-VASC is designed to store medical supplies related to vascular surgery and diseases of the vascular system.

The MS8140-VASC is preconfigured as follows: left column, top to bottom, 1ea catheter slide shelf, 1ea 1/2" shelf, 2ea 8" trays; right column, top to bottom, 3ea 2"trays, 5ea 4" trays, 3ea 8" trays.



Advantages and Benefits:

- Offers 35% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique "clip-on" locking design to withstand high loads.
- Unique design of plastic side panels allows modular trays to be extended and viewed at an angle without falling out of the cart. (pull-out stoppers required)
- Catheter slide shelf can be added featuring three cath slides, and five hooks per slide.

Standard Features:

- Catheter Slide Shelf featuring three slides and five hooks per slide
- One half inch shelf: Total of one shelf that is 23.50 x 15.75 x 1/2 inch
- Two inch modular tray: Total of three trays that are 23.50 x 15.75 x 2 inches
- Four inch modular tray: Total of five trays that are 23.50 x 15.75 x 4 inches
- Eight inch modular tray: Total of five trays that are 23.50 x 15.75 x 8 inches
- Tray dividers to organize width and depth of modular trays are sold separately
- Plastic side panels form a modular system consisting of 94 insertion opportunities (47 per column) to hold modular exchange trays and baskets.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with side of cart
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Sand, durable hard-baked hammer-beige powder coat textured finish
- Ships fully assembled with twelve year warranty

Assembled Cart Dimensions:

- 80.5" H x 42" W x 29.75" D
(204.47 cm H x 106.68 cm W x 75.565 cm D)
- Cart net weight: 247 pounds (112.04 kg)
- Approximate shipping weight: 295 pounds (133.81 kg)