



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Full Shell Medical Storage Cart Equipped with Five Shelves and Tambour Door #MS8140SH

The MS8140SH is a full shell medical storage cart equipped with five shelves and tambour door designed for quick access to medical instruments and supplies, featuring the most innovative modular tray system on the market. The MS8140SH can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- The MS8140SH has an open shell with five removable shelves. The open shell design offers 17% more storage capacity when compared with leading competitor's product.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- The shelves, with five dividers per shelf, are designed to improve visibility, accessibility and organization of sterile supply storage in the operating room.
- Plastic side panels allow for additional shelves to be mounted.
- Catheter slide shelf can be added featuring three cath slides, and five hooks per slide.

Standard Features:

- Five adjustable shelves with five adjustable dividers per shelf. Dividers can be mounted on ½ inch increments
- Plastic side panels form a modular system consisting of 47 insertion opportunities for mounting shelves
- Integrated push handles are flush with side of cart.
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Sand, durable hard-baked hammer-beige powder coat textured finish
- Ships fully assembled with twelve year warranty

Assembled Cart Dimensions:

- 80.5" H x 42" W x 29.75" D
(204.47 cm H x 106.68 cm W x 75.565 cm D)
- Cart net weight: 287 pounds (130.18 kg)
- Approximate shipping weight: 337 pounds (152.86 kg)