



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Preconfigured Triple Column Medical Storage Cart with Tambour Door #MS8160-CARDIO

The MS8160-CARDIO is a preconfigured triple column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8160-CARDIO is a preconfigured triple column medical storage cart designed to supplies for cardiovascular surgery to treat complications of congenital heart disease, rheumatic heart disease and atherosclerosis.

The MS8160-CARDIO is preconfigured as follows: left column, top to bottom, 3ea ½" shelves, 2ea 2" trays, 2ea 4" trays, 3ea 8" trays; middle column, top down, 4ea 2" trays, 3ea 4" trays, 4 ea 7" wired baskets; right column, top to bottom, 3ea ½" shelves, 2ea 2" trays, 2ea 4" trays, 4ea 8" trays.

Advantages and Benefits:

- Offers 35% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique "clip-on" locking design to withstand high loads.
- Unique design of plastic side panels allows modular trays to be extended and viewed at an angle without falling out of the cart. (pull-out stoppers required)
- Catheter slide shelf can be added featuring three cath slides, and five hooks per slide.

Standard Features:

- One half inch shelves: Total of six shelves that are 23.50 x 15.75 x ½ inch
- Two inch modular trays: Total of eight trays that are 23.50 x 15.75 x 2 inches
- Four inch modular trays: Total of seven trays that are 23.50 x 15.75 x 4 inches
- Seven inch wired baskets: Total of four baskets that are 23.50 x 15.75 x 7 inches
- Eight inch modular trays: Total of eight trays that is 23.50 x 15.75 x 8 inches
- Plastic side panels form a modular system consisting of 141 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Catheter Shelf Holders can be added to any column.
- One of the columns can be removed to increase catheter storage capacity
- Or remove column to install full-width shelves on the left or right of the cabinet. See product numbers MS8160SH-L or MS8160SH-R.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with side of cart.
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Sand, durable hard-baked hammer-beige powder coat textured finish
- Ships fully assembled with twelve year warranty



Assembled cart dimensions:

- 80.5" H x 59.625" W x 27.875" D
(204.47 cm H x 151.45 cm W x 70.803 cm D)
- Cart net weight: 375 pounds (170.10 kg)
- Approximate shipping weight: 425 pounds (192.78 kg)