



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Double Wide Full Shell Catheter Procedure Cart with Tambour Door Designed for use in a Cath Lab #MS8160CATH

The MS8160CATH Catheter Procedure Cart is a triple column cart with tambour door specifically designed for use in a cath lab, or an examination room in a hospital where diagnostic imaging equipment is used to support the catheterization procedure



Close Up Of Cath. Slide Shelf

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- The significant amount of catheter sizes and packaging techniques has made organizing and managing this critical area an extremely difficult task. Harloff offers a solution with a high density Catheter Storage System, which provides quick access and retrieval of catheters and related supplies.
- Equipped with three full extension catheter slides, each slide has five hooks, and each hook can hold up to 15 catheter packages. Each slide offers full extension, which brings the catheters into full view for quick access and visual inventory of remaining stock.
- The MS8160CATH stores long stent boxes, balloons, guidewires, and other supplies efficiently for easy access with our specially-engineered catheter and interior side panels.
- The maximum package length is 61 inches.

Standard Features:

- Catheter Slide Shelves, featuring nine cath slides and five hooks per slide, can be mounted in a lower position to create storage above for over-sized items.
- Cath slides consist of a pair of ball bearing drawer slides to withstand the pull force of removing catheters from hooks.
- Roll-up tambour door with key lock allows access to supplies during clinical hours
- Integrated push handles are flush with side of cart.
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Sand, durable hard-baked hammer-beige powder coat textured finish
- Ships fully assembled with twelve year warranty

Assembled Cart Dimensions:

- 80.5" H x 59.625" W x 27.875" D
(204.47 cm H x 151.45 cm W x 70.80 cm D)
- Cart net weight: 375 pounds (170.10 kg)
- Approximate shipping weight: 425 pounds (192.78 kg)