



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Triple Column Medical Storage Cart with Tambour Door #MS8160

The MS8160 is a triple column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market.

The MS8160 can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



Assembled cart dimensions:

- 80.5" H x 59.625" W x 27.875" D
(204.47 cm H x 151.45 cm W x 70.803 cm D)
- Cart net weight: 375 pounds (170.10 kg)
- Approximate shipping weight: 425 pounds (192.78 kg)

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Offers 35% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique "clip-on" locking design to withstand high loads.
- Unique design of plastic side panels allows modular trays to be extended and viewed at an angle without falling out of the cart. (pull-out stoppers required)
- Catheter slide shelf can be added featuring three cath slides, and five hooks per slide.

Standard Features:

- Plastic side panels form a modular system consisting of 141 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Catheter Shelf Holders can be added to any column.
- One of the columns can be removed to increase catheter storage capacity
- Or remove column to install full-width shelves on the left or right of the cabinet. See product numbers MS8160SH-L or MS8160SH-R.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with side of cart.
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Sand, durable hard-baked hammer-beige powder coat textured finish
- Ships fully assembled with twelve year warranty