

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Three Drawer Mini Line Infection Control Cart with Extended Stabilizing Bumper and 5" Wheels, with Key Lock #3243K

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints



Sixteen available color choices: *Antimicrobial Finish Available

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White

Specifications

Construction:

· Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters one brake, one directional
- Pre-threaded accessory mounting holes
- · Ball-bearing full extension drawer slides
- Removable plastic top
- · Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) 6-inch (15.2 cm)
- (2) 9-inch (22.9 cm)

Internal Drawer Dimensions:

• 12.6" W x 16.75" D (32cm x 42.5cm)

Assembled Cart Dimensions:

- 37" H x 19.5" W x 25" D (94.0 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 94 pounds (42.6 kg)
- Approximate shipping weight: 114 pounds (51.7 kg)