

## Mini Line Emergency Cart / Model: 3245B

Five Drawer with Breakaway Lock, Standard Package

Harloff's Mini Line of carts offers a narrow design with reduced footprints. Suitable for hospitals, clinics, and doctor offices with both space and budget constraints



### Seventeen Available Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Construction:

- Painted steel cabinet and drawers

### Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 24" (61 cm) of vertical drawer space
- Painted steel break-away lock bar with 100 break-away seals
- Stabilizing pontoon bumpers
- 5" (12.7 cm) ball-bearing caster set, all full swivel, two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

### Drawer Configuration:

- (3) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 9" (22.9 cm)

### Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5 cm)

### Overall Cart Dimensions:

- 35.75" H x 23.25" W x 22.2" D (90.8 cm H x 59.1 cm W x 56.4 cm D)
- Writing surface: 34.75" H x (88.3 cm)
- Cart net weight: 84 lbs (38.2 kg)
- Approximate shipping weight: 113 lbs (51.4 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (131.1 cm H x 73.7 cm W x 60.9 cm D)

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts