

Mini Line Isolation Cart/ Model: 3254K Four Drawer with Key Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- Stabilizing pontoon bumpers
- 5" ball-bearing full swivel casters one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) 3-inch (7.6 cm)
- (3) 9-inch (22.9 cm)

Internal Drawer Dimensions:

• 12.6" W x 16.75" D (32 cm x 42.5cm)

Overall Cart Dimensions:

- 43.75" H x 23.5" W x 22" D (110.5 cm H x 59.6 cm W x 55.8 cm D)
- Writing surface: 41.75" H x (106 cm)
- Cart net weight: 87 lbs (39.5 kg)
- Approximate shipping weight: 120 lbs (48.5 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (132 cm H x 73.6cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts