

Mini Line Anesthesia Cart/ Model: 3256K-ANS

Six Drawer with Key Lock, Specialty Accessory Package

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- Stabilizing pontoon bumpers
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-tier raised back rail system (MINI-RBR)
- (2) Aluminum side rails/push handles (68161)
- 6-bin tilt bin organizer (683805)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Overall Cart Dimensions:

- 58.75" H x 19.5" W x 25" D (149.2 cm H x 49.5 cm W x 63.5 cm D)
- Writing surface: 33.5" H x (81.2 cm)
- Cart net weight: 107 lbs (48.5 kg)
- Approximate shipping weight: 127 lbs (57.6 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (132 cm H x 73.6cm W x 60.9 cm D)