



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Five Drawer Mini24 Emergency Cart (Short Cabinet) with 3” Casters,
 Breakaway Lock
 #4145B**

Harloff’s Mini24 Line offers a 24” wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.



Color Shown is Red.
 Additional Colors Available.

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Breakaway lock system with package of 100 plastic seals
- 3” (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5” W x 16.75” D (47 cm x 42.5cm)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Assembled Cart Dimensions:

- 34.4” H x 24” W x 18”D
 (87.4 cm H x 61 cm W x 45.7 cm D)
- Cart net weight: 85 pounds (38.56 kg)

Shipping Dimensions:

- 49” H x 33” W x 26” D
 (124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 124 pounds (56.25 kg)