



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Five Drawer Mini24 Anesthesia Cart with Bumper and 5” Wheels,
 Key Lock
 #4255K**

Harloff’s Mini24 Line offers a 24” wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices



Color Shown is Light Blue.
 Additional Colors Available.

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 30 vertical inches (76.2 cm) of drawer space
- Key lock with two keys
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5” (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5” W x 16.75” D (47 cm x 42.5cm)

Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Shipping Dimensions:

- 57” H x 33” W x 26” D
 (144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 156 pounds (70.8 kg)

Assembled Cart Dimensions:

- 42.75” H x 29.25” W x 22” D
 (108.4 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 117 pounds (53.1kg)