

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

Classic Line Emergency Cart/ Model: 6331

Five Drawer with Breakaway Lock, Specialty Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

Sturdy, uni-body steel construction

Features:

- Stainless steel break-away lock bar
- Pack of 100 break-away seals
- 24" of vertical drawer space
- 5" (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, two braking
- Cabinet and drawer design allows for customized drawer configurations
- Painted steel top rail / push handle
- Antimicrobial plastic top, removable
- Stabilizing pontoon bumpers
- · Ball-bearing full extension drawer slides
- · Pre-drilled threaded holes for side mounted accessories
- · Cart is delivered mostly assembled
- Latex free
- · Durable, hard-baked powder coat finish

Accessories Included:

- Cardiac board and brackets (680406)
- Drop shelf (680300)
- Utility hooks (4 ea) (680402)
- IV pole (680401)
- Oxygen tank holder D/E size cylinders (AL680403)
- 6 compartment drawer trays (3 ea) (fit three per drawer), (40108)
- 7 compartment drawer tray (3 ea) (fit three per drawer), (40109)

Drawer Configuration:

- (3) 3-inch (7.6 cm)
- (1) 6-inch (15.2 cm)
- (1) 9-inch (22.5 cm)

Internal Drawer Dimensions:

• 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall Cart Dimensions:

- 39.5" H x 35" W x 22" D (100.3 cm H x 88.9 cm W x 55.9 cm D)
- Writing surface: 32.25" (89.5 cm)
- Height with IV Pole 79.5" (202 cm)
- Cart net weight: 131 lbs (59.4 kg)
- Approximate shipping weight: 170 lbs (69 kg)
- Approximate shipping dims: 47" H x 33" W x 26" D (119.4 cm H x 83.8 cm W x 66 cm D)