



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Nine Drawer, Pediatric Resuscitation Cart, Standard Package #6400PEC

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Tall cabinet with separate breakaway locks on each drawer
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cart is designed to be compatible with popular pediatric color coding systems
- Includes a pack of multi-colored labels for the end user to apply to identify drawers
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Anti-microbial plastic top, removable
- Painted steel top rail/ push handle
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled
- Twelve year warranty on manufactured parts

Drawer Configuration:

- (1) – 6" (15.2 cm)
- (8) – 3" (7.6 cm)

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 45.75" H x 32" W x 22" D
(116.2 cm H x 81.3 cm W x 55.9 cm D)
- Writing surface: 41.25" (104.8cm)
- Cart weight: 144 lbs (65.3 kg)
- Approx shipping weight: 185 pounds (83.9 kg)
- Approx. shipping dims: 57" H x 33" W x 26" D
(144.8cm H x 83.8cm W x 66cm D)

Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

