

Toll: (800) 433-4064 Ph: (719) 637-0300 medicalsales@harloff.com www.harloff.com

Classic Line, Six Drawer, Emergency Cart, Specialty Accessory Package



#### **Internal Drawer Dimensions:**

•23" W x 17" D (58.4 cm W x 43.2 cm D)

#### **Assembled Cart Dimensions:**

- 45.75" H x 35" W x 22" D (116.2 cm H x 89 cm W x 55.9 cm D)
- With IV pole: 86.25" (219 cm H)
- Cart weight: 150 pounds (68 kg)
- Approximate shipping weight: 189 pounds (85.7 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D (134 cm H x 83.8 cm W x 66 cm D)

# **Specifications**

#### **Construction:**

· Sturdy, uni-body steel construction

### Features:

- · Painted steel break-away lock bar
- · Pack of 100 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Painted steel top rail
- Stabilizing pontoon bumpers
- Five-inch (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- · Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- · Cart is delivered assembled
- Twelve year warranty on manufactured parts

# **Quick Ship:**

QuickShip applies to 6401Q red only, in quantities of 5 or less-Ships within 5 business days after receipt of order

#### **Accessories Included:**

- Cardiac board and brackets (680406)
- Drop shelf (680300)
- Utility hooks (4 ea) (680402)
- IV pole (680401)
- Oxygen tank holder D/E size cylinders (AL680403)
- 6 compartment drawer trays (3 ea) (fit three per drawer), (40108)
- 7 compartment drawer tray (3 ea) (fit three per drawer), (40109)

# **Drawer Configuration:**

- (4) 3" (7.6 cm) (1) 6" (15.2 cm)
- (1) 12" (30.5 cm)