



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 6 Drawer Workstation (Tall) with Key Lock, Deluxe Package #6456



**Color Shown is Hammertone Blue.  
Additional Colors Available.**

### Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

### Assembled Cart Dimensions:

- 42.5" H x 44" W x 22" D  
(108 cm H x 111.8 cm W x 55.9 cm D)
- With backrail: 57" H (144.8 cm H)
- Cart net weight: 160 lbs (72.5 kg)
- Approximate shipping weight: 201 lbs (91.2 kg)
- Approximate shipping dims: 53" H X 33" W X 26" D  
(134.6 cm H x 83.8 cm W x 66 cm D)

## Specifications

### Construction:

- Sturdy, uni-body All Steel Construction

### Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel  
- one tracking, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Stabilizing pontoon bumpers
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts
- Latex free

### Specialty Package Accessories Included:

- Raised back rail system with two horizontal rails (68163)
  - Four aluminum side rails / push handles (68161)
  - 6-Compartment tilt bin organizer (683805)
  - Tape and label dispenser (680415)
  - Catheter box, 18", with rail clip (680413)
  - Large utility basket (680407)
  - Plastic waste container, 3 gallon, with bracket and rail clip (684801, 683801)
  - Large sharps container and bracket with rail clip (684802, 683801)
  - 3-inch exchange trays (2) (681501)
  - 6-inch exchange trays (2) (681502)
  - 12-inch Adjustable drawer divider set (680503)
- Hospital Grade electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL specification change. Outlet conforming to UL1361A can be purchased as an accessory. (UL-OUTLETA)

### Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)