



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Three Drawer Infection Control Cart with Key Lock, Standard Package

#6510

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Aluminum side rail / push handle / Customer mount
- Stabilizing pontoon bumpers
- Cart is delivered fully assembled
- Twelve year warranty on manufactured parts
- Latex free
- Durable, hard-baked powder coat finish

Quick Ship:

QuickShip applies to 6510Q
Yellow only, in quantities of 5 or less-
Ships within 5 business days after receipt of order

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D
(108 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 112 pounds (50.8 kg)
- Approximate shipping weight: 148 pounds (67.2 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.6 cm H x 83.8 cm W x 66 cm D)



Color Shown is Yellow
Additional Colors Available

Isolation Cart Upgrade:

- Lemon drop yellow finish available (antimicrobial)