



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Classic Line, Three Drawer Infection Control Cart with Key Lock,
 Specialty Package
 #6511**



**Color Shown is Yellow
 Additional Colors Available**

Isolation Cart Upgrade:

- Lemon drop yellow finish available (antimicrobial)

Assembled Cart Dimensions:

- 45.75" H x 32" W x 22" D (116.2 cm H x 81.3 cm W x 55.9 cm D)
- With coat and hat rack: 63.5" H (161.3 cm H)
- Cart net weight: 114 pounds (63.6 kg)
- Approximate shipping weight: 150 pounds (68 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D (134.6 cm H x 83.8 cm W x 68 cm D)

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Removable plastic top
- Painted steel top rail / push handle (68160)
- Coat and hat rack (680412)
- Utility hooks (4 ea) (680402)
- Coat and hat rack / customer installed
- Cart is delivered mostly assembled
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

Quick Ship:

QuickShip applies to 6511Q
 Yellow only, in quantities of 5 or less-
 Ships within 5 business days after receipt of order

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)