

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

4 Drawer Workstation (Short) with Basic Electronic Pushbutton Lock, **Specialty Package** #7351E



Assembled Cart Dimensions:

- 36.5" H x 32" W x 22" D (92.7 cm H x 81.3 cm W x 55.9 cm D)
- With Back Rail: 51.75" H (131.5 cm H)
- Cart net weight: 112 lbs (50.8 kg)
- Approximate shipping weight: 152 lbs (69 kg)
- Approximate Shipping dims: 47" H x 33" W x 26" D
- (119.4 cm H x 83.8 cm W x 66 cm D)

Specifications

Construction:

• Sturdy, uni-body steel construction

Features:

- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 24" of vertical drawer space
- · Cabinet and drawer design allows for customized drawer configurations
- Four aluminum side rails / push handles
- · Raised back rail system with two horizontal rails
- · Ball-bearing full extension drawer slides
- · 5" ball-bearing casters, all full swivel one tracking, one with brake
- Removable plastic top
- · Pre-drilled threaded holes for side mounted accessories
- Stabilizing pontoon bumpers
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

Drawer Configuration:

- (2) 3" (7.6 cm) (1) 6" (15.2 cm)
- (1) 12" (22.9 cm)

Internal Drawer Dimensions:

• 23" W x 17" D (58.42 cm W x 43.18 cm D)