



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 6 Drawer Workstation (Tall) with Basic Electronic Pushbutton Lock, Specialty Package #7451E



**Color Shown is Beige.  
Additional Colors Available**

### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Raised back rail system with two horizontal rails
- Four aluminum side rails / push handles
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Stabilizing pontoon bumpers
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

#### Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

#### Assembled Cart Dimensions:

- 42.5" H x 32" W x 22" D  
(108 cm H x 81.3 cm W x 55.9 cm D)
- With Back Rail: 57" H (148 cm H)
- Cart net weight: 144 lbs (65.3 kg)
- Approximate shipping weight: 185 lbs (83.9 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D  
(134.6 cm H x 83.8 cm W x 66 cm D)