

## Aluminum Mini Line Emergency Cart/ Model: AL3256B

Six Drawers Aluminum with Breakaway Lock, Standard Package

Mini-Line Carts offer one of our smallest footprints. This cart has a 5" caster and provides more stability for carts that are moved frequently. Available in multiple heights, drawer configurations, and locking systems.



### Seventeen Available Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Construction:

- Painted aluminum cabinet and drawers

### Features:

- Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- Cabinet with 30" (76.2 cm) of vertical drawer space
- Painted steel break-away lock bar with 100 break-away seals
- Stabilizing pontoon bumpers
- 5" (12.7 cm) ball-bearing caster set, full swivel, two with brake
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

### Drawer Configuration:

- (3) – 3" (7.6 cm)
- (2) – 6" (15.2 cm)
- (1) – 9" (22.9 cm)

### Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5 cm)

### Overall Cart Dimensions:

- 42.75" H x 23.25" W x 22" D (108.6 cm H x 59.1 cm W x 55.9 cm D)
- Writing surface: 41.75" H x (106 cm)
- Cart net weight: 76 lbs (34.5 kg)
- Approximate shipping weight: 104 lbs (47.3 kg)
- Approximate shipping dims: 52" H x 29" W x 24" D (132 cm H x 73 cm W x 58.4 cm D)

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts