



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Aluminum Universal Line Tall 6 Drawer Anesthesia Cart, with Basic Electronic Pushbutton Lock and Anesthesia Accessory Package #AL809E6-ANS

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.



Color Shown
is Beige.
Additional Colors
Available

Drawer Configuration:

- (4) – 3.25" (8.3 cm)
- (1) – 6.5" (16.5 cm)
- (1) – 9.75" (24.8 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 66.5" H x 32.5" W x 24.75" D
(169 cm H x 83 cm W x 63 cm D)
With Back Rail System
- Cart weight: 117 pounds (53.2 kg)
- Approximate shipping weight: 147 pounds (67 kg)

Specifications

Construction:

- Aluminum construction

Features:

- Internal pullout shelf on left side
- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- 29.25" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- 5" (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered wheels-on and mostly assembled
- Pre-drilled hole locations for standard accessories – many are pre-threaded for ease of installation
- Latex free
- Durable, hard-baked powder coat finish available in 16 colors
- Twelve year warranty on manufactured parts.

Package Accessories Include:

- 4 Aluminum Side Rails
- Three Tier Back Rail System
- Hospital Grade Electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL Specification change. Outlets conforming to UL1363A can be purchased as an accessory item. (Part #UL-OUTLETA)

OPTIONAL UPGRADE: ANS2 PACKAGE

AL809E6-ANS2

Includes all ANS package items plus:

- Six compartment tilt bin
- Label holder
- Catheter box
- Three gallon sharps container
- Three gallon waste container
- Three inch full drawer tray with dividers
- Two (2) each, six inch exchange trays
- Nine inch drawer divider set