

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 \* Ph: (719)637-0300 Fax:(719)597-8273 \* www.harloff.com

# Universal Line Aluminum Super 6 Drawer Cart with Breakaway Lock and **Deluxe Emergency Package** #AL810B6-EMG2

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with

flexible configuration options.



### **Construction:**

Aluminum construction

#### Features:

- · Internal Pullout shelf on right side side
- Break-away lock bar and 100 numbered break-away seals
- 32.5" of vertical drawer space
- · Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- · Removable plastic top
- Extended caster footprint for stability with protective bumper
- 5" (12.7 cm) ball-bearing premium casters, all full swivelone tracking, one braking
- · Ball-bearing full extension drawer slides
- · Cart is delivered wheels-on and mostly assembled
- · Pre-drilled hole locations for standard accessories many are pre-threaded for ease of installation
- Latex free
- Durable, hard-baked powder coat finish available in 17 colors
- · Twelve year manufacturer's warranty

# Package Accessories Include:

- · Adjustable Back Shelf with Hanger Bar
- Suction Shelf
- · Cardiac board and brackets
- Utility hooks (2 ea)
- IV pole
- Aluminum oxygen tank holder D/E size cylinders
- Full Drawer Trays with adjustable divider set, 68530-P1 (1 ea)
- (4) 3" Exchange Trays with dividers
- (2) ea. 6" Exchange Trays with dividers
- Hospital Grade Electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL Specification change. Outlets conforming to UL1363A can be purchased as an accessory item. (Part #UL-OUTLETA)

# **Drawer Configuration:**

- (3) 3.25" (8.3 cm) (2) 6.5" (16.5 cm) (1) 9.75" (24.8 cm)

# **Internal Drawer Dimensions:**

• 23" W x 17" D (58.4 cm W x 43.2 cm D)

#### **Assembled Cart Dimensions:**

- 52" H x 45.5" W x 25" D
- (132 cm H x 115.6 cm W x 63.5 cm D)
- IV pole extends to 88.5" (225 cm) H
- Cart weight: 145 pounds (66 kg)
- Approximate shipping weight: 188 pounds (85.5 kg)



Addtional Colors Available