



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

MR-Conditional, Seven Drawer Emergency Cart, Breakaway Lock, Specialty Package #MR7B-EMG



Color Shown is Red. Additional Colors Available
Color Choice References Drawers,
Cabinet Remains White

Emergency Specialty Package Includes MRI-Compatible Accessories:

- Stainless steel IV pole (MR-IV2)
- Plastic cardiac board and painted aluminum mounting brackets (MR-CARDBRD)
- Painted aluminum oxygen tank holder – D/E size cylinders (MR-ALO2)
- Aluminum utility hooks (4 ea.) (MR-UHOOKS4)
- Plastic full drawer trays with adjustable dividers (2 ea.) (68530-P1)

Construction:

- Non-magnetic or “weakly magnetic” materials including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware

Features:

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional”
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Low-ferrous 5” sealed ball-bearing casters; full swivel – one with brake, one directional
- Painted aluminum breakaway locking system
- Breakaway lock seals, numbered, yellow, pkg-100
- Full wrap-around vinyl bumper
- Stainless steel pull-out shelf
- Replaceable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Twelve year warranty on manufactured parts
- Hard-baked powder coat finish

Drawer Configuration:

- (1) - 4” (10.2 cm)
- (4) - 3.25” (8.3 cm)
- (1) - 6.5” (16.5 cm)
- (1) - 9.75” (24.8 cm)

Internal Drawer Dimensions:

- 23.9” W x 17” D (60.8 cm x 43.2 cm)

Assembled Cart Dimensions

- 45.75” H x 30.5” W x 23.75” D
(116.2 cm H x 77.5 cm W x 60.3 cm D)
- Cart net weight: 118 pounds (53.5 kg)
- Approx. shipping weight: 135 pounds (61.2 kg)