



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 3/4 Height Preconfigured Double Wide Catheter Procedure Cart with Tambour Door

### #MS6140-CATH3

The MS6140-CATH3 is a ¾ height full shell, double wide Catheter Procedure Cart with tambour door, specifically designed for use in a cath lab to support the catheterization procedure. At five feet tall, the cart is 110 pounds lighter and over a foot shorter than the MS8140CATH

The MS6140-CATH3 comes preconfigured as follows: left column, top to bottom 1 catheter slide shelf with 3 slides; right column, top to bottom, 3ea 4" trays, 3ea 8" trays



### Assembled Cart Dimensions:

- 60.5" H x 42" W x 27.9" D  
(153.7 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 240 pounds (108 kg)
- Approximate shipping weight: 370 pounds (167 kg)

### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- The significant amount of catheter sizes and packaging techniques has made organizing and managing this critical area an extremely difficult task. Harloff offers a solution with a high density Catheter Storage System, which provides quick access and retrieval of catheters and related supplies.
- Equipped with adjustable Cath Slide Shelf, consisting of seven full extension catheter slides, each slide has five hooks, and each hook can hold up to 15 catheter packages.
- Each Cath Slide features a pair of ball bearing slides to withstand the pull force of removing catheters from hooks; full extension brings the catheters into full view for quick access and visual inventory of remaining stock.
- Cath Slide Shelf can be mounted in top position or 13 inches lower for storage above.
- Stores long stent boxes, balloons, guidewires, and other supplies efficiently for easy access with our specially-engineered catheter.
- The maximum package length is 55 inches.

### Standard Features:

- Catheter Slide Shelf, featuring three cath slides and five hooks per slide, can be mounted in a lower position to create storage above for over-sized items. Cath Slide Shelf consists of a pair of ball bearing drawer slides for extra strength and durability
- Four inch modular trays: Total of three trays that are 23.50 x 15.75 x 4 inches
- Eight inch modular trays: Total of three trays that are 23.50 x 15.75 x 8 inches
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty