



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Single Column Medical Storage Cart with Tambour Door #MS8128

The MS8128 is a single column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market.

The MS8128 can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using any 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular panels to be viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains eight hooks.

### Standard Features:

- Plastic side panels form a modular system consisting of 47 insertion opportunities to hold modular exchange trays and baskets, which are sold separately.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Standard color is Sand. Also available in 16 other colors
- Ships fully assembled with twelve year warranty

### Assembled Cart Dimensions:

- 80.5" H x 24.4" W x 29.75" D  
(204.5 cm H x 69.9 cm W x 75.5 cm D)
- Cart net weight: 190 pounds (86.2 kg)
- Approximate shipping weight: 308 pounds (139.7 kg)