



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Preconfigured Double Column Medical Storage Cart with Tambour Door MS8140-ARTHRO

The MS8140-ARTHRO is a preconfigured double column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8140-ARTHRO Cart is designed to support arthroscopic surgery by storing supplies required for procedures performed using an arthroscope, a type of endoscope that is inserted into the joint through a small incision.

The MS8140-ARTHRO is preconfigured as follows: left column, top to bottom, 3ea ½ " trays, 2ea 4" trays, 1ea 8" tray, 1ea 7" wired basket; right column, top to bottom, 1ea 1/2" tray, 1ea 2" tray, 5ea 7" wired baskets.



### Assembled Cart Dimensions:

- 80.5" H x 42" W x 27.9" D  
( 204.5 cm H x 106.7 cm W x 75.8 cm D)
- Cart net weight: 345 pounds (156.5 kg)
- Approximate shipping weight: 475 pounds (215.5 kg)

### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be extended and viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains eight hooks.

### Standard Features:

- One half inch trays: Total of four trays that are 23.50 x 15.75 x 1/2 inches
- Two inch modular trays: Total of one tray that is 23.50 x 15.75 x 2 inches
- Four inch modular trays: Total of two trays that are 23.50 x 15.75 x 4 inches
- Seven inch wired baskets: Total of six baskets that are 23.50 x 15.75 x 7 inches
- Eight inch modular trays: Total of one tray that is 23.50 x 15.75 x 8 inches
- Tray dividers to organize width and depth of modular trays are sold separately
- Plastic side panels form a modular system consisting of 94 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty