



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Preconfigured Double Column Medical Storage Cart with Tambour Door #MS8140-A

The MS8140-A is a preconfigured double column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8140-A can be used in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.

The MS8140-A is preconfigured as follows: right column, top to bottom, 2ea 2" trays and 6ea 8" trays; left column, top to bottom, 2ea 2" trays and 11ea 4" trays



Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using any 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains five hooks.

Standard Features:

- Two inch modular trays: Total of four trays that are 23.50 x 15.75 x 2 inches
- Four inch modular trays: Total of eleven trays that are 23.50 x 15.75 x 4 inches
- Eight inch modular trays: Total of six trays that are 23.50 x 15.75 x 8 inches
- Plastic side panels form a modular system consisting of 94 insertion opportunities (47 per column) to hold modular exchange trays and baskets.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty

Assembled Cart Dimensions:

- 80.5" H x 42" W x 27.9" D
(204.5 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 330 pounds (149.7 kg)
- Approximate shipping weight: 460 pounds (208.7 kg)