



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Preconfigured Double Column Medical Storage Cart with Tambour Door #MS8140-B

The MS8140-B is a preconfigured double column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8140-B can be used in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.

The MS8140-B is preconfigured as follows: left column, top to bottom, 4ea 4" wired baskets, 5ea 7" wired baskets; right column, top to bottom, 4ea 4" wired baskets, 5 ea 7" wired baskets.



### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- Offers 30% more tray or basket storage capacity by volume when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains five hooks.

### Standard Features:

- Four inch wired baskets: Total of eight baskets that are 23.50 x 15.75 x 4 inches
- Seven inch wired baskets: Total of ten baskets that are 23.50 x 15.75 x 7 inches
- Plastic side panels form a modular system consisting of 94 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty

### Assembled Cart Dimensions:

- 80.5" H x 42" W x 27.9" D  
( 204.47 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 318 pounds (144.2 kg)
- Approximate shipping weight: 448 pounds (203.2 kg)