



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Preconfigured Triple Column Medical Storage Cart with Tambour Door #MS8160-CYSTO

The MS8160-CYSTO is a triple column cart with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8160-CYSTO is designed to supplies for surgical procedures related to the urinary bladder and the urinary tract.

The MS8160-CYSTO is preconfigured as follows: left column, top to bottom, 1ea 4" tray, 6ea 8" trays; middle column, top down, 13ea 2" trays, 1ea 8" tray; right column, top to bottom, 1ea catheter slide, 1ea 8" tray.



Assembled cart dimensions:

- 80.5" H x 59.6" W x 27.9" D
(204.5 cm H x 151.4 cm W x 70.8 cm D)
- Cart net weight: 400 pounds (181 kg)
- Approximate shipping weight: 590 pounds (267.2 kg)

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be extended and viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains eight hooks.

Standard Features:

- Two inch modular trays: Total of thirteen trays that are 23.50 x 15.75 x 2 inches
- Four inch modular trays: Total of one tray that is 23.50 x 15.75 x 4 inches
- Eight inch modular trays: Total of seven trays that is 23.50 x 15.75 x 8 inches
- Plastic side panels form a modular system consisting of 141 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Catheter Shelf Holders can be added to any column.
- One of the columns can be removed to increase catheter storage capacity
- Or remove column to install full-width shelves on the left or right of the cabinet. See product numbers MS8160SH-L or MS8160SH-R.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side, positioned vertically for better control
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty