



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Standard Line 294 Capacity Unit Dose Medication Box Cart, Standard Package #SL420BOX

### Specifications



#### Construction:

- Painted steel cabinet and aluminum drawers

#### Features:

- Dual column medication cart, single width narcotics drawer
- Capacity of 420 unit-dose boxes, 42 per drawer
- BEST® brand lock on cabinet; Illinois lock on narcotics drawer
- Painted steel pull-out side shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Offered in specialty package (**#SL420BOXSP**) which includes: BEST® lock on narcotics drawer, #68640-P1 drawer tray with adjustable dividers, #680551 eye drop tray, AL680435 aluminum waste container and #684810 locking sharps container

#### Thirteen available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY) Hammertone Green (-HTGN) and Sand (-SD).

#### Assembled Cart Dimensions

- 45.75" H x 43" W x 23.5" D  
(116.2 cm x 109.2 cm x 59.7 cm)
- Writing Surface: 44" (111.8 cm)
- Cart weight: 230 lbs (104.6 kg)
- Approx shipping weight: 275 lbs (125 kg)

#### Drawer configuration:

##### Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 9 ea – 3.25" med box drawers with internal divider set, 42 box capacity per drawer

##### Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer