



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## V-Series, Six Drawer Procedure Cart with Key Lock #V30-6K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.

### Specifications



#### Construction:

- Steel construction

#### Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

#### Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

#### Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D  
 (101.6 cm H x 74.9 cm W x 55.9 cm D)
- Cart weight: 118 pounds (53.5 kg)
- Approximate shipping weight: 141 pounds (63.9 kg)

#### Drawer Configuration:

- (4) – 3" (7.6 cm)
  - (1) – 6" (15.2 cm)
  - (1) – 12" (30.5 cm)
- Integrated Pull Handle on all drawers

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black