



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Value Line Medication Cart, Basic Electronic Pushbutton Lock with Key Lock Override(MEDICINE ON TIME® compatible)

VLT150EMOT

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- Capacity of up to 150 Medicine On Time® cards
- Internal locking narcotics box made for M.O.T. packaging in bottom 9” drawer
- Pull out side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5” ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified “Latex Free”
- Five year warranty



Color Shown is Beige



Drawer Configuration:

- 1 - ea 3” multipurpose storage drawer
- 2 - ea 9.75” drawers with 1 adjustable divider
- 1 - ea 9.75” drawer with locking narcotics compartment made for M.O.T. packaging.

Internal Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 46.75” H x 29.5” W x 22.5” D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25” (109.9 cm)
- Cart weight: 140 pounds (63.7 kg)
- Approx. shipping weight: 175 pounds (79.6 kg)

Thirteen available color choices:

(may be viewed at www.harloff.com)

- Sand, Gloss Beige, Burgundy, Light Blue, Yellow, Navy, Light Gray, Mauve, Red, Teal, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green