



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## Value Line Medication Cart with Key Lock # VLT216BOX

Harloff’s value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.



Color Shown is Light Gray



### Thirteen available color choices:

(may be viewed at [www.harloff.com](http://www.harloff.com))

- Beige, Burgundy, Light Blue, Dark Blue, Light Gray, Mauve, Red, Teal, Yellow, Hammertone Blue, Hammertone Gray, Hammertone Green and Sand

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Tall cabinet with key locking system with 4 keys
- Internal locking narcotics box in bottom 9” drawer
- Capacity of up to 216 unit dose boxes or 60 sections for strip packaging
- Pull out Side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5” ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified “Latex Free”
- Twelve year manufacturer’s warranty

#### Drawer Configuration:

- 1 - ea 3” multipurpose storage drawer
- 6 - ea 3” drawers with internal unit dose dividers
- 1 - ea 9” drawer with 1 divider and locking narcotics box

#### Internal Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

#### Assembled Cart Dimensions:

- 46.75” H x 29.5” W x 22.5” D (121.3 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 43.25” (109.9 cm)
- Cart weight: 170 pounds (77.3 kg)
- Approx. shipping weight: 210 pounds (95.5 kg)