



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Wooden Laminate In Wall Medication Cabinet, Single Door, Single Lock #WL2786-CM

Add the look of wood to functional storage for an attractive in-room storage solution. Purpose designed to meet industry standards for security and function. Cabinet is designed for in wall stud mounting. Slam lock ensures security with one hand open and close

Specifications

Construction:

- Steel cabinet
- Door is MDF panel covered with a special heat-pressed thermal foil laminate with a rich wood appearance

Features:

- In-room medication cabinet, capacity for 20 punch cards/liquid storage
- Wood laminate door in Dark Cherry or Solara Oak
- Door opens to create an area for medication preparation
- Slam lock provides increased assurance that unit is locked
- Multiple cabinets may be keyed differently or keyed alike –Must be specified on order
- In wall mount with hardware included
- Medication cabinet is delivered ready to use
- Twelve year warranty against manufacturer defect
- Durable hard-baked powder coat finish paint on box – Color is gloss beige.



Dark
 Cherry
 (DC)

WL2786-DC



Solara
 Oak
 (SO)

WL2786-SO

Color Options:

The WL2786-DC is shown in Dark Cherry finish also available in Solara Oak

Additional laminate color options available with minimum cabinet purchase and extended lead time.

Overall Dimensions:

- 13.25" H x 15.5" W x 6.75" D
 (33.7 cm H x 39.4 cm W x 17.1 cm D)
- Net weight: 14.5 pounds (6.6 kg)
- Approx. shipping weight: 17 pounds (7.7 kg)