

Wood Vinyl 6 Drawer Emergency Cart, Standard Package #WV6400-CM

Hospitals, Long Term Care facilities, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's WV line of vinyl-coated aluminum carts. Designed and built to support culture change initiatives, this cart gives the buyer the best of both worlds, combining the aesthetics of wood with the long-term functionality of a Harloff metal cart.



Construction:

- Vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate

Features:

- This cart is designed to look like a piece of furniture
- Light weight vinyl coated aluminum
- Break-away lock bar with 100 break-away seals
- 5" (12.7cm) ball-bearing premium casters, all full swivel, 1 braking, 1 directional
- Ball-bearing full extension drawer slides
- Pull out shelf
- Cart is washable with soap and water
- Twelve year manufacturer's warranty

Drawer Configuration:

- 4 ea - 3.25" drawers
- 1 ea - 6.5" drawer
- 1 ea - 12" drawer

Assembled Cart Dimensions:

- 44.5" H x 30" W x 22.5" D
- Writing surface: 44"
- Cart weight: 115 pounds
- Approx. shipping weight: 150 pounds

Color Choices:

- The WV6400-CM is shown above in Cherry Mahogany finish with a brown granite finish on the top. Also available in Buffalo Oak.
- Other laminate colors available but require minimum purchase and extended lead times