



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## OptimAL Line, Aluminum 540 Punch Card Medication Cart #AL540PC



### Specifications

#### Construction:

- Painted aluminum cabinet and drawers

#### Features:

- Single column medication cart with 540 punch card capacity, single width internal narcotics box
- BEST® brand lock on cabinet. Illinois lock on narcotics drawer
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Painted steel pull-out shelf
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#AL540PCSP) which includes #68640-P1 drawer tray with adjustable dividers, AL680436 aluminum waste container with lid and #684810 locking sharps container

#### Thirteen available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY) Hammertone Green (-HTGN) and Sand (-SD).

#### Assembled Cart Dimensions

- 45.75" H x 36" W x 23.5" D  
(116.2 cm x 91.4 cm x 59.7 cm)
- Writing Surface: 44" (111.8)
- Cart weight: 137 lbs (62.2 kg)
- Approx shipping weight: 197 lbs (89.6 kg)

#### Internal Drawer Dimensions:

- 18" D x 27.5" W

#### Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea - 9.75" punch card drawer with dividers and locking narcotics box
- 2 ea – 9.75" punch card drawers with dividers