



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

**Standard Line 540 Punch Card Medication Cart, Standard Package  
 #SL540PC**



**Specifications**

**Construction:**

- Painted steel cabinet and aluminum drawers

**Features:**

- Single column medication cart with 540 punch card capacity, single width internal narcotics box
- BEST® brand lock on cabinet. Illinois lock on narcotics drawer
- Painted steel pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Offered in specialty package (**#SL540PCSP**) which includes: #68640-P1 drawer tray with adjustable dividers, AL680435 aluminum waste container and #684810 locking sharps container

**Thirteen available color choices:**

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY), Hammertone Green (-HTGN) and Sand (-SD).

**Assembled Cart Dimensions**

- 45.75" H x 34" W x 23.5" D (116.2 cm x 86.4 cm x 59.7 cm)
- Writing Surface: 44" (111.8)
- Cart weight: 142 lbs (72.8 kg)
- Approx shipping weight: 185 lbs (93.2 kg)

**Internal Drawer Dimensions:**

- 18" D x 27.5" W

**Drawer Configuration:**

- 1 ea – 4" multipurpose storage drawer
- 1 ea - 9.75" punch card drawer with dividers and locking narcotics box
- 2 ea – 9.75" punch card drawers with dividers