



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Value Line Medication Cart with Key Lock # VLT360PC

Specifications



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system
- Taller drawer height for increased capacity
- Capacity of up to 360 punch card carts
- Internal locking narcotics box in bottom 9" drawer
- Painted steel pull out side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension ball-bearing drawer slides
- Cart is delivered fully assembled
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish
- Certified "Latex Free"
- Five year warranty

Available in "SP" Specialty Package which includes locking sharps container, waste container and top drawer tray with adjustable dividers. #VLT360PCSP

Thirteen available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-HTDB), Hammertone Grey (-HTGY), Hammertone Green (-HTGN) and Sand (-SD).

Assembled Cart Dimensions

- 46.5" H x 29.5" W x 22.5" D
(118 cm H x 75 cm W x 57.2 cm D)
- Writing surface: 42.5" (108 cm)
- Cart weight: 140 pounds (63.7 kg)
- Approx. shipping weight: 175 pounds (79.6 kg)

Drawer Dimensions:

- 23" W x 17" D (58.5 cm W x 43.2 cm D)

Drawer Configuration:

- 1 x 3" multipurpose storage drawer.
- 2 x 9.75" drawers with two adjustable dividers.
- 1 x 9.75" drawer with 1 divider and locking narcotics box.