

FULL LINE CATALOG



HARLOFF
INNOVATING STORAGE FOR CARE


ABOUT US



Harloff has built a reputation for creating a personalized experience that's designed to meet your specific needs. We specialize in innovative designs that are standardized to the medical profession, as well as building custom designed products to meet your unique requirements.

We have an in-house engineering staff that will work with you to design specialized products, as well as outside "in the field" consultants who bring various levels of knowledge and expertise to all of our products. Our products are proudly "American made" by our staff here in Colorado Springs.

Harloff has taken great care in cultivating a national network of local expert representatives that provide a personalized experience and fast accurate service for all of our customers. We value your business and take an attitude of ensuring that all your needs are met in a friendly, professional, and timely manner.



Our diverse production capabilities include working with expansive resources and materials, which include stainless steel, painted steel, aluminum, wood laminate and wood vinyl.

Harloff has been creating innovative carts and storage to meet your needs since 1951. We take pride in our team and look forward to providing you with an enriching customer experience.

Harloff is indeed Innovating Storage for Care, helping you deliver quality care for your patients.

Using This Catalog:

We manufacture the finest quality carts and storage systems to suit every aspect of healthcare storage, security and mobility. Our sales and distribution partners are carefully chosen to work with you on the solutions that specifically meet your needs. Our catalog represents only the base of our product lines and is meant to help you work with us to match the exact solution that fits your needs.

Core Products –

Anesthesia Carts & Workstations.....	3
Crash Carts.....	4
Infection Control & PPE.....	5
Medication & Pharmacy Carts.....	6
Treatment & Procedure Carts.....	7
Cast & Splinting Carts.....	8
MRI Solutions.....	9
Specialty Products.....	10

Storage Products –

Casework.....	11
Laminate Wood Storage.....	12
Scope Storage.....	13
Catheter Storage.....	14
Supply Storage.....	15
Lean Storage Systems.....	16

Value Solutions, Accessories and Quick Ship –

Mini & Value Carts.....	17
Quick Ship Options.....	18
Cart Accessories.....	19
Security Options.....	20
Standard Features and Color Chart.....	21
Ordering Information.....	22

Set Up Where You Need It

We manufacture the highest quality anesthesia storage carts in a variety of sizes that are ideal for use in surgical environments and improve operating room utilization through increased efficiency that improve the patients experience. Some of our Anesthesia Carts are made of “non-magnetic” or “weakly magnetic” materials including aluminum, stainless steel, plastic, and aluminum and brass hardware, and are safe to use in MRI suites.

Features & Benefits:

- Predrilled holes for mounting accessories
- Protective pontoon bumpers
- Internal pull-out shelf
- Toprail/push handle
- Multiple locking options
- Ball-bearing, full extension drawer slides
- 5” premium ball-bearing casters, full swivel – one brake, one directional



Universal Tall
Anesthesia Cart
#AL809E6-ANS



IV and Epidural
Start Cart
#3123



Classic Line
6-Drawer Cart
#6450



Mini24
Anesthesia
Cart
#4256E

Accessories



Divider Set
680503



Divider Tray
68530-PI



Waste Bin
684801



Tilt Bins
683803



Back Rail
68163

Designed For Emergency Environments

Harloff has 4 standard product families to meet your budget, space and accessory needs for emergency care. We are also well known for our ability to design custom carts to meet your specific needs. Multiple security systems are available, along with multiple drawer configurations, and a full line of cart accessories for unique facility requirements and specialized applications.



Tall
6-Drawer
#6401

Features & Benefits:

- Sturdy unibody construction
- Stabilizing pontoon bumpers
- 5" ball-bearing premium casters, all full swivel, one tracking and one braking
- Removable plastic top
- Ball-bearing full extension drawer slides
- Full line of accessories available
- 3 Quick ship models available



Short
5-Drawer
#3145B



Lightweight
Aluminum Construction
5-Drawer
#AL808B5



Short
4-Drawer
#V24-4B

Accessories



O² Tank Holder
AL680403



IV Pole
680401



Cardiac Board
680406



Drawer Trays
40108 & 40109



Divider Tray
68530-P1

Designed to Help Prevent Infection

Intended to help prevent the spread of infectious illness and disease within a healthcare setting. Our carts are designed to store all supplies required by protocol to reduce the risk of pathogen transmission between healthcare workers and patients. Helps lower hospital costs and liability by reducing the number of hospital acquired infections.

Features & Benefits:

- Silver ion paint coating with antimicrobial properties available with the #ISO6550 model upon request, in lemon drop yellow only
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel, one directional, one with toe-locking brake
- Cabinet and drawer design allows for customized drawer configurations
- Aluminum side rail / push handle available as an accessory



Open Well
Isolation
Station
#ISO6550



4-Drawer
#6520



Mini
4-Drawer
#3154K



Tall
4-Drawer
#V30-4K

Wall mount or over the door

Saves floor space and allows storage of personal protective equipment in areas with space limitations.



ISO2800

ISO2835

ISO2825

Secure Medication Storage

Harloff medication carts are manufactured to provide hospital pharmacies with highly secure medication administration in a healthcare environment. Our carts are manufactured with durable, high capacity drawers, which improve medication management and ensure patient safety while reducing costs.

Features & Benefits:

- Compatible with most accessories, with pre-threaded accessory mounting holes
- Locking narcotics options available for all carts
- Full extension ball-bearing drawer slides
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Pull-out shelf included as standard
- Divider options for all medications packaging



Medication Cart
#VLT24BIN5



Punch Card Cart
#SL600PC



Medication Bin Cassette Cart
#AL16CS-5



Wood Look Medicine Cabinet
#WL2780RECTRM-DC

Accessories



Laptop Arm
680327



Waste Bin
684801



Sharps Container
684810



Waste Container
AL680436



Cup / Straw Holder
AL2360

Easy Access To Supplies

Treatment carts are the perfect solution for storing all the supplies needed in one cart. Our carts help keep numerous small supplies organized for efficient and effective application that improve patient outcomes and experience.

Features & Benefits:

- Sturdy uni-body steel or aluminum construction
- The #7450E comes with an Electronic Pushbutton Lock
- Stabilizing pontoon bumpers
- 5" ball-bearing premium casters, all full swivel – one brake, one tracking
- Ball-bearing full extension drawer slides
- Keyless entry options



6-Drawer
#7450E



Compact
5-Drawer
#V24-5K

6-Drawer Mini24
Space Saver
#4256K



Accessories



Turns any Room into a Casting Area

Harloff makes it easy to store and transport instruments and supplies for the treatment of strains, sprains and fractures. These storage carts are designed for use in emergency departments, urgent care centers, ambulatory surgery centers and hospital based orthopedic surgery services. Designed in conjunction with orthopedic technicians to maximize functionality.

Features & Benefits:

- Stainless steel bucket holder with vinyl mat
- Aluminum towel bar
- 24 inch stainless steel ruler on top of cart
- Full wrap-around vinyl bumper on separate stabilizing frame
- Internal drawer dividers
- 5" ball-bearing premium casters, all full swivel – one tracking, one with brake
- Ball-bearing, full extension drawer slides



8-Drawer Stainless Steel Cart with Deluxe Accessory Package and Top Compartment
#6025-TC



8-Drawer with Deluxe Accessory Package
#6035



Standard Package 8-Drawer Casting Cart
#6030

Accessories



Shelf Insert
602002



Top Compartment
602106



Tourniquet Pole
602007



Scissors Cable
602004



Folding Shelf
602001

Made for Imaging

Harloff is the only USA based manufacturer to offer MR-Conditional Safe carts in both 1.5 and 3.0 Tesla magnetic resonance environments. Our carts are built primarily of non-ferrous aluminum, plastic and stainless steel for safe daily use in highly magnetic environments. We also offer packages for Anesthesia and Emergency requirements, multiple drawer configurations, and extra MR conditional accessories.

Features and Benefits:

- Specially sourced, low ferrous, ball-bearing full extension drawer slides
- BEST® brand key lock with two keys
- 5" MR conditional full swivel casters
- White outer cabinet, with 17 color choices for drawers
- MRI decal that identifies the cart as belonging to MRI suites



6-Drawer Standard Cabinet with Specialty Package
#MR6K-MAN



6-Drawer Narrow Body Cabinet
#MRN6K

6-Drawer Standard Cabinet
#MR6B



Accessories



Tilt Bins
MR-5TILT/4TILT



Sharps Container
MR-3GSHARPS



Cardiac Board
MR-CARDBRD



Waste Bins
MR-8QWASTE



O² Bracket
MR-ALO2



IV Pole
MR-IV2

All of Harloff's Specialty Carts feature the same superior construction and specialized features as our other carts, and are designed for specialized treatments and uses.



#IOL924

INTRAOCCULAR LENSES

Allows you to perform more procedures each day by storing lenses in one mobile unit. Compatible with multiple brands and sizes of IOL boxes, including Alcon®, Bausch and Lomb®, Rayner®, and Abbott® Medical Optics. Provides secure storage for expensive equipment and supplies.

MALIGNANT HYPERTHERMIA

Designed with a high performance Follett® REF2 medical grade refrigerator that ensures temperature control for cold isotonic saline for IV infusion and/or gastric peritoneal or rectal irrigation that helps you ensure patient safety during emergency situations. Our MH cart meets the standards for JACHO Best Practices.



#MH5200B



#6400PEC

PEDIATRIC RESUSCITATION

Be ready to care for your smallest patients. This cart is designed to be compatible with popular pediatric color coding systems, and separate breakaway locks for each drawer. Accessory package available that includes cardiac board, drop shelf, utility hooks, IV pole, and oxygen tank holder.

NARCOTIC BOXES

Our narcotic storage boxes offer maximum security - constructed with 18 and 20 gauge steel, high security pick-resistant locks, and full length steel pinned door hinges to prevent pin removal, with predrilled holes in the back for mounting to a solid surface. The #2810A and #2820A come with an anti-microbial finish.



#2703

Double Door/Double Lock
- Heavy Duty



#2810A

#2820A

Both Single & Double Door
Locks - in a Variety of Sizes



#27AVD20

Double Door/Double Lock
- with Audio Visual Alarm



2729E

Single Door/Single
Electronic Push-Button
Lock

Configurable Style, Durability, and Aesthetically Pleasing

Harloff Laminate medical casework provides the optimum in modular and flexible static storage options for lean supply management in the OR, Lab, storage rooms, and procedure rooms. Designed to look like built-in mill work that can be combined in a custom fashion to accommodate any workspace, while featuring high density, flexible storage interior, adjustable trays, bins and baskets. Save space, organize, and increase efficiency for improved accessibility of supplies resulting in better quality of care.

Features and Benefits:

- Modular, reconfigurable cabinets
- Custom grade MDF panels covered with high pressure laminate exterior and white low pressure laminate interior
- Flexible side panels for modular shelves or tray system
- 270 degree opening 3-knuckle hinge



Combine for a customized look



Wall Cabinet
#MSWCS3035-GD



Base Cabinet
#MSFCB3435-SD
(Counter top not included)



Single Vertical Cabinet
#MSFCB8118-GD

Casework Laminate Colors



Accessories

WORKS WITH ALL

MEDSTOR^{MAX}

ACCESSORY PRODUCTS

Add Wood For Warmth

Harloff provides an attractive alternative for facilities looking to soften their institutional look. Our medication furniture is designed to meet industry standards for security and function.

Wood provides an attractive alternative to traditional medication and treatment carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of your hospital, clinic or long term care facility.



Drawer & Cabinet Treatment Cart
#WV6200-CM

Features & Benefits:

- Washable with soap and water
- Adjustable drawer divider options for medications, liquids, eye drops, and other items
- Ball-bearing full extension drawer slides
- All carts are moveable (models with skirting covers 4" ball-bearing casters)



4-Drawer Punch-card Medication Cart
#WL480PC-DC



4-Drawer Narrow Punch-card Medication Cart
#WL240PCN-SO



6-Drawer Emergency Cart
#WV6400-CM



Shown with Optional Side Cabinet

Storage Laminate Colors



Dark Cherry and Solara Oak are Standard Laminate Colors for Carts and Cabinets Ordered in any Quantity, OTHER COLORS REQUIRE MINIMUM PURCHASE.

MEDSTORMAX[®]

Cabinets for Storing Scopes

All of our state-of-the-art stationary scope storage cabinets meet AAMI ST91 standards. We offer an optional state-of-the-art HEPA filter with a removable drip tray, and your choice of Tambour or glass doors. Our interiors are engineered to protect and safely store scopes while providing easy access.

Features & Benefits:

- Constructed of sturdy 18 gauge steel
- Our large cabinet features 2 drying fans 120VAC (220 with international plug adapter #INTLPLUGADPTR)
- Optional HEPA[®] filter captures up to 99.97% of particles, 0.3 microns, exceeding AORN scope storage recommendations; with removable drip tray
- Large cabinets hold up to 16 scopes in a secure and stationary location
- Small cabinets hold up to 8 Bronchoscopes
- The large scope cabinet comes with choice of key locking double doors with tempered glass, or roll-up tambour door



Scope Storage with Tambour Door with Drying Package
#SC16DP

Tall Scope Storage with Glass Doors with Drying Package
#SC16DRDP

Short Scope Storage with Tambour Door
(HEPA filter not included)
#SCBR-8



SCOPE STORAGE

MEDSTOR^{MAX}

Catheter Accessibility

Harloff catheter storage cabinets come either mobile or fixed, and offer high density secure storage for valuable medical supplies used in Cardiology. Our storage cabinets are designed for quick access to medical supplies, featuring the most innovative modular tray system on the market, and a choice of two different door styles for added security and ease of access. Our storage carts are superior in quality, actual storage space, and cost to competitive offerings.



Single Column Catheter Rack Storage with Glass Doors
#MS8128DRCATH

Double Column Catheter Rack Storage Cabinet with Glass Doors
#MS8142DRCATH



Features & Benefits:

- Single, double, and triple column carts can be configured in numerous ways to meet your special needs
- Cath Lab in a cart: designed especially for Cath Labs and IR departments
- Holds all sizes of catheters on full extension holders
- All-inclusive carts that keeps all the necessary supplies together in one place
- Cath slides feature ball-bearing slides to withstand the pull force of removing catheters from hooks
- Catheter shelf can be added featuring 3 robust slide assemblies, each containing 5 hooks

Double Column Split Storage & Catheter Racks w/Tambour Door
#MS6140CATH4



Accessories



RACK LABELS



CATH RACK 3



CATH RACK 7

CATHETER STORAGE

MEDSTORMAX[®]

Maximum Flexibility

Harloff's Medstormax medical supply storage carts are designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. Our carts are either mobile or stationary, and offer a choice of two different door styles for security and ease of access, providing high density and secure storage for valuable medical supplies.

Features & Benefits:

- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be extended and viewed at an angle without falling out of the cart
- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4", and 8" modular trays; or 4" and 7" wire baskets
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect" to withstand high loads
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase costs



Triple Column Storage
with Tambour Door
#MS8160-ORTHO

Double Column Storage
with Glass Doors
#MS8140DR-A



3/4 Height Storage Cart
#MS6141DR-SH

Accessories



WIRE BASKETS



PULL-OUT SHELF



PLASTIC TRAYS

SUPPLY STORAGE

OPEN FRAME

Lean Management of Inventory

Open storage solutions provides a clean, easy to see, and a systematic replenishing system that makes managing inventory efficient, accurate, and timely. It improves functionality of an organization and reduces labor and waste.




Features & Benefits:

- Easy to see and determine stocked inventory levels, color coded by category and basket labeling available
- Works very well as a two-bin Kanban system
- Min/Max inventory system
- Standard PAR Level inventory management process
- Our storage equipment achieves the highest density of storage capacity per square foot

WIRE SYSTEMS

Easy Visibility, Durable and Clean

Time-proven wire systems can be configured a variety of ways to accommodate specific inventory needs while helping implement lean management systems by improving functionality and reducing waste. These units are durable and easy to clean.



Features & Benefits:

- Frames are strong, durable, easy to keep clean, and last a long time
- Easy visibility for maintaining and ordering supplies
- All inventory is easy to access, improves processes and outcomes, reduces costs, and increases efficiency
- Shelves can be configured in numerous ways to meet your specific stocking needs
- Shelves are modular and flexible



Economical and Saves Space

Our Mini-line and Value carts are cost-effective and our space-saver design can be used in environments with limited space to safely secure supplies.



#3254E



#3154K



#3156E



#3254B



#4156E-ANS



#4255K



#4256B



#4736K



#VLT12BIN5



#VLT360PC



#V30-6B



#V24-4EMG



#V30-2K



#V30-3K



#V24-3K

Full line of accessories available - Call today for more information and pricing

Fast Availability

Though we don't stock pre-made items, we offer a variety of products that can be shipped quickly to meet urgent needs. All of our Quick-ship models can be constructed and shipped within five business days from receiving the order. Because of this fast turnaround time, they're only available in the colors shown, and any other standard colors may take longer to ship.



#6301Q



#6400Q



#6401Q

#VLT360PCQ



#SL600PCHBG



#2701BG



#2702BG



#2703BG



#2730BG



#2820AQ



#2810AQ



Quick Ship boxes available in Beige only. Models #2820AQ and #2810AQ finished in anti-microbial paint.

Full line of accessories available - Call today for more information and pricing

Add Accessories to Customize your Cart

Accessories allow you to customize your cart to meet specific needs or to improve utilization in a variety of medical environments.



O² Bracket
#AL680403



IV Pole
#680401



Cardiac Board
#680406



Protected Outlet Strip
#680404A



Waste Container w/wire
#684801



8 qt. Waste Bin
#680408



Sharps Container/bracket
#684802



Drawer Dividers
#680503



Glove Box Holder
#AL2362



Cup & Straw Holder
#AL2359



Notebook Holder
#AL2366



Waste Container
#AL680435



Defib Shelf
#680323



Locking Sharps Holder
#684810



Tilt Bin
#683803



3" Exchange Tray
#681501



Adjustable Dividers
#680574



3" Drawer Divider Tray
#68530P1



7/6 Compartment Trays
#40109 & #40108



6" Exchange Tray
#681502

*Multiple Tray Sizes Available. Contact Harloff for More Information

Multiple Locking Options

Harloff understands the importance of safely securing medical supplies in a public, hospital environment, and offers a variety of high-quality standard locks, closure options, and key lock overrides. We can also use custom mechanisms for special security needs upon request.



STANDARD KEY LOCK



BASIC ELECTRONIC PUSH-BUTTON LOCK WITH KEY LOCK OVERRIDE



BREAKAWAY LOCK



WIRELESS ELECTRONIC PUSH-BUTTON LOCK WITH KEY LOCK OVERRIDE



EKC FLUSH KEYPAD WITH KEY LOCK OVERRIDE



DOUBLE DOOR/DOUBLE LOCK WITH ALARM & LIGHT



SINGLE DOOR/DOUBLE LOCK





SINGLE DOOR/SINGLE LOCK

SECURITY OPTIONS

Our Standard Features:



Standard Colors:



Textured Colors:



Customer Service hours:

Monday - Friday

10 am - 7 pm Eastern

9 am - 6 pm Central

8 am - 5 pm Mountain

7 am - 4 pm Pacific



Placing your order

- Orders can be received by fax, mail, email, or sent EDI with established partners.
- No phone orders will be taken.
- Orders for painted product must specify a color to be considered complete and entered for production. Please see Colors on page 21.
- Order acknowledgments will be sent within 24 hours to confirm details of your order, i.e. products ordered, pricing and estimated ship date.
- All changes to purchase orders must be confirmed by hard copy.

Terms of sale

- Net 30 days for established customers.
- Visa, MasterCard, American Express, ACH and wire transfers accepted.
- C.O.D. shipments are not available.

Product warranties

- 12-year limited warranty on parts manufactured by Harloff.
- 1-year limited warranty on parts and hardware purchased by Harloff.
- Please contact customer service for additional details.

On-time delivery

Harloff takes pride in having a superior on-time delivery rate, and quick turn-around times for all your product orders. One of our primary goals is to ship orders in a timely manner to meet promised delivery dates. Our 95% on-time delivery target is not only met, but often exceeded. Harloff employs the principles of lean manufacturing to improve efficiency and on-time delivery, while constantly striving to get quality products out the door with the shortest lead times possible.

In a rapidly changing world, Harloff offers the creativity and flexibility needed to accommodate your specific needs in a consistent and predictable manner.

HARLOFF

INNOVATING STORAGE FOR CARE

650 Ford Street
Colorado Springs, CO 80915

719-637-0300

800-433-4064

sales@harloff.com

www.harloff.com

 Contract Holder

