



MEDICATION CARTS **AL300BOX**

UNIT-DOSE MEDICATION BOX CART

300 Box Capacity, 2-Column, Standard Line

Features:

- Painted aluminum cabinet and drawers
- Dual column medication cart with 300 unit dose medication box capacity, single width narcotics drawer
- BEST® brand lock on cabinet with four keys
- Illinois lock on narcotics drawer with four keys
- Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Heat sealed powder coat finish
- Dividers are available to accommodate different medication packaging
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Color shown is Hammertone Blue. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Standard Line 300 Capacity Unit Dose Medication Box Cart

AL210BOX medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.

Left Column Drawer Configuration

- 1 ea - 4" multipurpose storage drawer
- 9 ea - 3.25" drawers with internal dividers

Right Column Drawer Configuration

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 9.75" locking narcotics drawer
- 2 ea - 9.75" multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D (116.2 cm H x 91.4 cm W x 59.7 cm D)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 167 pounds (76 kg)
- Approx. shipping weight: 207 pounds (94.1 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm x 119.4 cm x 68.6 cm)



View similar products online at

<https://www.harloff.com/product-category/multi-dose-medication-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300